


Using LabVIEW to Automate the ALEPH High Intensity Titanium Sapphire Laser



Sunny Wong, Reed Hollinger, Jorge Rocca
sw758@cornell.edu
Colorado State University REU Summer 2020



Project Overview

- 1 Examples of experiments that could benefit from automation
2. A basic data acquisition loop
3. Interfacing with Newport Power Meters for energy measurements
4. Interfacing with stepper motors to control the camera and target system
5. Interfacing with the laser firing system.


Benefits of Automation

- Current shot rates: one shot per 1-5 minutes
- The 0.85-petawatt Ti:Sa Laser can run at 3.3 hertz
- With enough data points through intensity scans, diagnostics behave as a function of energy

A	B	C	D	E	F	G	H	I	J	K
	Time 1	Shot Count	Energy 1 (mJ)	Energy 2 (mJ)	Target	Comments	Post-Comments	Prepulse Energy	Prepulse Delay	Prepulse Duration
	7/21/2020 6:42:43 PM	1	0.393	19.99558	25 um Ni	Starting experimental scan on Ni foil observing K shell		~10 mJ	0.5 ns	60 fs
	7/21/2020 6:49:57 PM	2	0.361	19.71816	25 um Ni	Starting experimental scan on Ni foil observing K shell	K-a still > He-a	~10 mJ	0.25 ns	60 fs
	7/21/2020 6:54:43 PM	3	0.35	19.07796	25 um Ni	Starting experimental scan on Ni foil observing K shell		~10 mJ	0.125 ns	60 fs
	7/21/2020 7:00:27 PM	4	0.391	19.16332	25 um Ni	Starting experimental scan on Ni foil observing K shell		~10 mJ	62.5 ps	60 fs
	7/21/2020 7:06:55 PM	5	0.388	19.16332	25 um Ni	Starting experimental scan on Ni foil observing K shell		~10 mJ	30 ps	60 fs
	7/21/2020 7:11:33 PM	6	0.441	18.6725	25 um Ni	Starting experimental scan on Ni foil observing K shell		~10 mJ	single pulse no prepulse	60 fs
	7/21/2020 7:20:37 PM	7	0.368	18.3524	25 um Ni	Starting experimental scan on Ni foil observing K shell		~10 mJ	15 ps	60 fs
	7/21/2020 7:26:46 PM	8	0.37	18.43776	25 um Ni	Starting experimental scan on Ni foil observing K shell		~10 mJ	trying for 0 ps	60 fs
	7/22/2020 12:31:03 PM	9	0.218	10.3499	CH 4 micron	Espec/Proton Wedge/ Bubble Detectors/ NTOF (should be 10 mJ)		10 mJ	1 ns prepulse	60 fs
	7/22/2020 12:34:15 PM	10	0	11.54494	CH 4 micron	test energy shot		10 mJ	1 ns prepulse	60 fs
	7/22/2020 12:44:00 PM	11	0	8.77074	CH 4 micron	test energy shot		10 mJ	1 ns prepulse	60 fs
	7/22/2020 12:58:17 PM	12	0.406	16.30376	CH 4 micron	Shooting again with 1 ns delay	Shot 12 prepulse may have been 10 mJ	10 mJ	1 ns prepulse	60 fs
	7/22/2020 1:01:44 PM	13	0.404	18.43776	CH 4 micron	blocked prepulse		10 mJ	no prepulse	60 fs
	7/22/2020 1:07:56 PM	14	0.383	17.8189	CH 4 micron	1 ns prepulse again (may have been misaligned last shot)		10 mJ	1 ns	60 fs

Example: Neutron Generation


- High intensity laser irradiates an array of deuterated polyethylene (CD_2) nanowires to drive a D-D fusion reaction
- D-D fusion reaction generates a neutron burst
- The maximum number of neutrons per shot was 3.6×10^6 with 60 femtosecond laser pulses at 1.64 J
- Limited to shot/min, but high repetition rates are desirable for fusion science and high energy density physics




Curtis, A., et al. *Nature Communication* **9**, 1077 (2018).

Example: X-Ray Tomography

- Using X-Ray sources (created by irradiating a tantalum foil with high intensity laser) to generate tomographic images of electronics
- The GIF to the left consists of 72 shots (taken at about a shot per minute)
- An automated laser system running at 3.3 Hz can complete the same scan in about 20 seconds.




Data Acquisition Control Loop



Newport 1919-R Power Meters

1. Program receives measurements from both power meters




2. Plots measurements and monitors

Fundamental And Second Harmonic Energy

- Uses bar graphs and arrays to plot information in real time

3. Logs Data as a CSV file

- Allows user to add sub-headers, add comments, delete shots, and to change file-path



Graphical Interface

File Edit View Project Operate Tools Window Help

Save Data?

Shot Time

7/29/2020

5:28:00 PM

Last Logged Shot Count
85.000000

File Path to Save Data

C:\Users\Laser Lab\Documents\BackupVIs\logbook20200716.csv

Enter FilePath

Note: The file will not save properly if it is open while the procedure runs.

Target for Data

CD2 foam 150mg/cc 4um foam Type F

Enter Target

Comment Box for Data

Au coating

Enter Comment

Note: Writing a comment that contains commas will cause the string to be split apart and stored in separate columns.

Post-Comment

no electrons on shot likely out of focus

Enter Post-Comment

Note: Post comments will only be entered once per press, and the VI will write the string to the last recorded row in the CSV file.

Additional Column Headers Already in File (Indicator)

--	--	--	--

Additional Column Headers (Input)

--	--	--

Additional Comments

10 mJ single shot 60 fs

Enter Text

Enter Text

STOP

Shot Count To Delete
60

Delete Shot

Current Energy Reading 1
0.559

2nd Harmonic Side Window

Energy (mJ)

Plot 0

Multiplier for 2nd harmonic side window
1000


Current Energy Reading 2
0.12804

Fundamental Energy


Energy (mJ)

Plot 0

Multiplier for fundamental
213400




Code



Maneuvering the Camera System and Target Mount

1. LabVIEW program sends TTL signals to relay switches via NI CompactDAQ Digital I/O board
2. Relay Switches are wired in parallel to the physical stepper motor system built by CSU
3. Program receives and displays response signal to confirm that the motor controls are moving



Firing the Laser

1. Check if cameras are down
2. If the above condition is satisfied, the laser will fire.

In the future, an interlocking camera system will help protect equipment.

More functionality will be added later.



Summary

The LabView program...

1. Displays and logs data from the two power meters
2. Can maneuver the camera and target system
3. Can fire the laser after checking if the camera is down

Email: Sw758@cornell.edu